

F5 Medical Bed

MLH-01 · Technical Datasheet

Full-featured Saudi-made medical bed with electric patient control, comprehensive nurse-side control panel, central brake system and rechargeable backup battery.

Key Features

- Patient control: head and knees, up/down
- Nurse control: lockouts, cardiac chair, Trendelenburg, auto-contour, CPR control
- Side rails with integrated controls
- Central braking and bumpers
- Rechargeable backup battery
- Integrated foot extension, head and foot panels, drainage bag holder



Specifications

Mattress	Pressure-reducing, 14 cm thick
Power	Mains + rechargeable backup battery
CPR	Nurse-side CPR control
Optional	Foam mattress, IV pole, night light, overbed tables (MLA-23/24)