


Multi-Functional PRIME care Bed

MLP-01 · Technical Datasheet

Multi-functional electric in-patient bed with stainless steel IV pole, high-endurance silent castors and detachable ABS head and foot boards, supplied with bumpers.

Key Features

- Multi-function electric positioning
- Stainless steel IV pole
- High-endurance silent castors with strong brakes
- Removable ABS head and foot panels
- Supplied with protective bumpers



Specifications

Length	200 cm
Width	90 cm
Height	81 cm
Optional	Backup battery, CPR system, radiolucent backrest, under-bed lighting, bed extension, central braking