

Santron Pro ICU Bed

SANTRON PRO · Technical Datasheet

The flagship of the MEDcare Santron family. Advanced ICU configuration with integrated weighing scale, digital nurse control panel, bed-exit alarm and enhanced monitoring for demanding ICU and long-term care. Manufactured in Saudi Arabia.

Key Features

- Integrated weighing scale
- Nurse control panel with display
- Footboard with display screen
- Bed-exit alarm with low-position indicator
- Under-bed and side lighting
- All standard Santron positioning and safety functions



Specifications

Length	215 cm
Width	95 cm
Height	45–75 cm
Weighing	Integrated patient scale
Controls	Digital nurse panel + side patient panel
Alarm	Bed-exit alarm with low-position indicator
Optional	Orthopedic arm