

Santron Hospital Bed

SANTRON · Technical Datasheet

HOSPicare's flagship electric ICU bed. Advanced clinical positioning (cardiac chair, Trendelenburg), integrated X-ray cassette holder, central braking and under-bed lighting. Hospital-grade, manufactured in Saudi Arabia.

Key Features

- Integrated X-ray cassette holder
- Side patient control panel and standard footboard control
- Clinical positioning modes: flat examination, cardiac chair, bed entry/exit, Trendelenburg and reverse Trendelenburg
- Central braking with directional lock
- Under-bed lighting and side bumpers with integrated IV pole
- Hygienic removable sleep surface, hospital-grade construction



Specifications

Length	215 cm
Width	95 cm
Height	45–75 cm
Positioning modes	Flat, cardiac chair, Trendelenburg, reverse Trendelenburg
Brakes	Central braking with directional lock
X-ray	Integrated cassette holder
Optional	Orthopedic arm